

The tradition  
of universal therapy.

**Qs**  
Magnetotherapy



**ASV**  
A S A L A S E R  
V E T E R I N A R Y  
Research and Therapeutic Solutions

# ASA Magnetotherapy for veterinary applications

The functioning of ASA Magnetotherapy devices, based on the scientific principles of ELF (Extremely Low Frequency 1–100 Hz) pulsed magnetic fields, is characterized by low frequency and low intensity.

- ▶ **Easy Qs Vet** is ideal for the specific treatment of localized areas thanks to the Flexa Vet applicators. Easy Qs Vet offers pre-set treatment programs, giving the operator the possibility to customize the emission parameters based on the type of animal, pathology and clinical phase.
- ▶ **PMT Qs Vet** inherits the features of Easy Qs Vet and broadens its applications: thanks to independent channels it allows the use of portable solenoids to treat the patients.



Discover more about  
Magnetotherapy  
Qs Vet:

## Advantages

- ▷ Action even on deep tissues
  - ▷ Well-tolerated
  - ▷ Non-invasive
  - ▷ Painless
  - ▷ Direct action on the whole body
  - ▷ Can be used as stand-alone therapy or in combination with other therapies
- ▽ **Easy Qs Vet**  
The portable solution
  - ▽ **PMT Qs Vet**  
The most complete and customizable version



# Examples of Application

## ▶ OSTEOARTICULAR LEVEL

- ▶ Fractures and delays in consolidation
- ▶ Osteoarthritis
- ▶ Arthritis

## ▶ VASCULAR & TISSUE REPAIR LEVEL

- ▶ Wounds and Ulcers
- ▶ Edemas

## ▶ NEUROMUSCULAR LEVEL

- ▶ Peripheral nerves lesions
- ▶ Muscle tears
- ▶ Contusions
- ▶ Sprains



# BioCote® antimicrobial silver ion technology

**CERTIFIED TECHNOLOGY THAT REDUCES THE LEVEL OF BACTERIA, MOULDS AND FUNGI ON SURFACES BY UP TO 99.99%.**



The BioCote® silver ion technology of the Flexa and Solenoid cover:

- ▶ is resistant to body fluids (blood and urine) and is easily cleaned with common detergents;
- ▶ offers 24-hour, continuous and integral applicators protection;
- ▶ retains its antimicrobial characteristics over time, without wearing out or dissolving;
- ▶ is hard-wearing and comfortable to the touch, providing a comfortable place for the animals to lie during therapy;
- ▶ the cover also features zips that allow 2 applicators to be joined together to form a larger mat.





**CORPORATE HEADQUARTERS / REGISTERED OFFICE**

Via Galileo Galilei, 23 / 36057 Arcugnano (VI) - Italy  
T +39 0444 28 92 00 / F +39 0444 28 90 80

[asavet@asalaser.com](mailto:asavet@asalaser.com)

**RESEARCH DIVISION / BRANCH**

Joint Laboratory Department of Experimental and Clinical  
Biomedical Sciences University of Florence  
Viale G. Pieraccini, 6 / 50139 Florence - Italy

[asacampus@asalaser.com](mailto:asacampus@asalaser.com)

[asaveterinary.com](http://asaveterinary.com)

**Copyright © ASA srl** - All rights reserved. Complete or partial copying, printing and distribution of the information present in this document is prohibited, unless written consent from the owner is granted. This content is intended for healthcare professionals and does not substitute a professional medical advice.

02/2025 - EN

